

---

# **In-Vivo-Imaging-Pipeline**

***Release 0.0.0***

**Darik A. O'Neil**

**Jan 03, 2023**



**CONTENTS:**

<b>1</b>	<b>Management</b>	<b>1</b>
1.1	Organization module . . . . .	1
1.2	UserInterfaces module . . . . .	10
<b>2</b>	<b>Indices and tables</b>	<b>11</b>
	<b>Python Module Index</b>	<b>13</b>



## MANAGEMENT

### 1.1 Organization module

**class** Organization.**BehavioralExperiment**(*Meta: Tuple[str, str], ExperimentName: str*)

Bases: *Experiment*

Data Class for a generic day of a behavioral task

#### Required Inputs

*Meta* : Passed meta from experimental hierarchy (directory, mouse\_id)

*ExperimentName* : Title of ExperimentName

#### Properties

*mouse\_id* : Identifies which mouse this data belongs to

*instance\_data* : Identifies when this behavioral ExperimentName was created

#### Attributes

*data* : Pandas dataframe of synced data

*folder\_dictionary* : A dictionary of relevant folders for this behavioral ExperimentName

*modifications* : List of modifications made to this behavioral ExperimentName

*meta* : bruker metadata

*multi\_index*: Pandas multi-index of behavioral components

*state\_index* : look-up table / index relating states to integers

*trial\_parameters* : behavioral parameters

#### Methods

*add\_image\_sampling\_folder* : Generates a folder for containing imaging data of a specific sampling rate

*load\_data* : Loads all data

*record\_mod* : Records a modification made to the behavioral ExperimentName (Date & Time)

*update\_folder\_dictionary* : This function re-indexes all folders in the folder dictionary

**load\_data**(*ImagingParameters: Optional[Union[dict, list[dict]]] = None, \*args: Optional[Tuple[str, str]], \*\*kwargs*) → Self

Loads all data

### Parameters

- **ImagingParameters** (*Optional[dict]*) – Parameters for some imaging dataset or list of datasets (e.g., for two different sampling rates)
- **args** (*Tuple[str, str]*) – Optionally pass Sync Key to synchronize bruker recordings
- **kwargs** – passed to internal functions taking kwargs

### Return type

Any

**class** Organization.Data(*Path: str*)

Bases: object

This is a class for managing a folder of unorganized data files

### Required Inputs

*Path* : path to folder

### Self Methods

*find\_matching\_files* : Finds all matching files

*reindex* : Function that indexed the files within folder again

*find\_all\_ext* : Finds all files with specific extension

### Properties

*instance\_data* : Data created

*path* : path to folder

*files* : List of files in folder

**property files:** List[str]

**find\_all\_ext**(*Ext: str*) → Optional[List[str]]

Finds all files with specific extension

### Parameters

**Ext** (*str*) – Filename extension

### Returns

List of files

### Return type

List[str]

**find\_matching\_files**(*Filename: str, Folder: Optional[str] = None*) → Optional[Tuple[str]]

Finds all matching files

### Parameters

- **Filename** (*str*) – Filename or ID to search for
- **Folder** (*Any*) – Specify folder filename in

### Returns

Matching file/s

**Return type**

Any

**property folders: dict**

Dictionary of folders in path

**Return type**

dict

**property instance\_date: str**

Date Created

**Return type**

str

**property path: str**

Path to folder

**Return type**

str

**reindex()** → Self

Function that indexes the files within folder again

**class** Organization.**Experiment**(*Meta: Tuple[str, str], ExperimentName: str*)

Bases: object

Class for a generic experiment

**Required Inputs***Meta* : Passed meta from mouse (directory, mouse\_id)*ExperimentName* : Title of experiment**Properties***instance\_data* : Identifies when this behavioral experiment was created*mouse\_id* : Identifies which mouse this data belongs to**Attributes***data* : a pandas dataframe containing synchronized data*folder\_dictionary* : A dictionary of relevant folders for this experiment*meta* : bruker metadata*modifications* : List of modifications made to this experiment**Public Methods***add\_image\_sampling\_folder* : Generates a folder for containing imaging data of a specific sampling rate*load\_data* : Loads all data*record\_mod* : Records a modification made to the experiment (Date & Time)*update\_folder\_dictionary* : This function re-indexes all folders in the folder dictionary**property ExperimentName\_id: str**

**add\_image\_sampling\_folder**(*SamplingRate: int*) → Self

Generates a folder for containing imaging data of a specific sampling rate

**Parameters**

**SamplingRate** (*int*) – Sampling Rate of Dataset in Hz

**Return type**

Any

**copy\_raw\_imaging\_data**() → Self

This function copies raw imaging data to the appropriate folder

**Return type**

Any

**property instance\_date: str**

Date created

**Return type**

str

**load\_data**(*ImagingParameters: Optional[Union[dict, list[dict]]] = None, \*args: Optional[Tuple[str, str]], \*\*kwargs*) → Self

Loads all data

**Parameters**

- **ImagingParameters** (*Optional[dict]*) – Parameters for some imaging dataset or list of datasets (e.g., for two different sampling rates)
- **args** (*Tuple[str, str]*) – Optionally pass Sync Key to synchronize bruker recordings
- **kwargs** – passed to internal functions taking kwargs

**Return type**

Any

**Return type**

Any

**property mouse\_id: str**

ID of mouse

**Return type**

str

**record\_mod**() → Self

Records a modification made to the behavioral ExperimentName (Date & Time)

**Return type**

Any

**update\_folder\_dictionary**() → Self

This function re-indexes all folders in the folder dictionary

**Return type**

Any

**class** Organization.**Figures**(*Path: str*)

Bases: *Data*

A class for storing figures, inherits collected data folder



**view\_figure**(*Name: str*) → plt.Figure

Function identifies and views a figure based on supplied name

**Parameters**

**Name** (*str*) – Name of figure (can be partial)

**Returns**

the plotted figure

**Return type**

Any

**class** Organization.**Images**(*Path: str*)

Bases: *Data*

Class specifically for folders containing raw images, inherits collected data folder

**property channels**

**property file\_format**

**property frames**

**property height**

**property imaging\_files**

**property meta\_files**

**property num\_imaging\_files**

**property num\_meta\_files**

**property planes**

**reorganize\_bruker\_files**() → None

This function extracts out the meta files and saves in a new directory

**Return type**

None

**property width**

**class** Organization.**ImagingAnalysis**(*Path: str*)

Bases: *Data*

Class specifically for imaging analysis folders, inherits collected data folder

**Self Methods**

*load\_fissa\_exports* : loads fissa exported files

*load\_cascade\_exports* : loads cascade exported files

*load\_suite2p* : loads suite2p exported files

*export\_registration\_to\_denoised* : moves registration to new folder for namespace compatibility when skipping denoising step

*clean\_up\_motion\_correction* : This function removes the reg\_tif folder and registered.bin generated during motion correction.

*clean\_up\_compilation* : This function removes the compiled tif files

*add\_notes* : Function adds notes

**add\_notes**(*Step: str, KeyOrDict: Union[str, dict], Notes: Optional[Any] = None*) → Self

Function adds notes indicating steps

**Parameters**

- **Step** – Step of Analysis
- **Step** – str
- **KeyOrDict** (*Union[str, dict]*) – Either a Key or a dictionary containing multiple key-value (note) pairs
- **Notes** (*Optional[Any]*) – If using key, then notes is the paired value

**Return type**

Any

**clean\_up\_compilation**() → Self

This function removes the compiled tif files generated inside CompiledImagingData (You can avoid the creation of these in the first place by changing suite2p parameters)

**Return type**

Any

**clean\_up\_motion\_correction**() → Self

**This function removes the reg\_tif folder and registered.bin generated during motion correction.**  
(You can avoid the creation of these in the first place by changing suite2p parameters)

**Return type**

Any

**property current\_ExperimentName: str**

ExperimentName of Analysis

**Return type**

str

**default\_folders**()

**export\_registration\_to\_denoised**()

moves registration to new folder for namespace compatibility

**Returns**

**load\_cascade\_exports**() → Tuple[ndarray, ndarray, ndarray, dict]

This function loads the Spike Times, Spike Prob, Discrete Approximation and ProcessedInferences files exported from Cascade

**Returns**

SpikeTimes, SpikeProb, DiscreteApproximation, Processed Inferences

**Return type**

tuple[Any, Any, Any, dict]

**load\_fissa\_exports**() → Tuple[dict, dict, dict]

This function loads the prepared and separated files exported from Fissa

**Returns**

Prepared, Separated, ProcessedTraces

**Return type**

tuple[dict, dict, dict]

**load\_suite2p**(\*args: str)**class** Organization.**Mouse**(\*\*kwargs)

Bases: object

Class for Organizing &amp; Managing Experimental Data Across Sessions

**Keyword Arguments***Logfile* : Path to existing log file (str, default None)*Mouse* : Mouse ID (str, default None)*Condition* : Experimental Condition (str, default None)*Directory* : Directory for hierarchy (str, default None)*Study* : Study (str, default None)*StudyMouse* : Study ID (str, default None)**Properties***mouse\_id* : ID of Mouse*log\_file* : Log Filename Path*experimental\_condition* : Experiment condition of the mouse*instance\_data* : Date when this experimental hierarchy was created**Attributes***directory* : Experimental Hierarchy Directory*experiments* : Names of included experiments*study* : Study*study\_mouse* : ID of mouse in study*modifications* : modifications made to this file**Public Class Methods***load* : Function that loads the entire mouse**Public Methods***create* : This function creates the directory/logs/organization.json if it doesn't exist*check\_log* : Checks Log Status*create\_log\_file* : Creates log file*pass\_meta* : Passes directory/mouse id*record\_mod* : Record modification of experiment*record\_experiments\_mod* : Record modification of experiments*save* : Saves mouse to organization.json*start\_log* : Starts Log**Private Class Methods***\_generate\_analysis\_subdirectory* : Generate Analysis

`_generate_analysis_technique_subdirectory` : Generate Analysis Technique  
`_generate_behavior_subdirectory` : Generate Behavioral Folder  
`_generate_directory_structure` : Generates the Directory Structure (The structured folders where data stored)  
`_generate_experiment_folders` : Generate Behavioral ExperimentName Folder  
`_generate_histology_directory` : Generates Histology Folder  
`_generate_imaging_subdirectory` : Generate Imaging Folder  
`_generate_roi_matching_index_directory` : Generate ROI Matching Folder

**check\_log()** → Self

Checks log status

**Return type**

Any

**create()** → Self

This function generates the directory hierarchy in one step

**Return type**

Any

**create\_experiment**(*ExperimentName: str, Type: Optional[str, object] = 'Experiment', \*\*kwargs*) → Self

Generates an experiment ExperimentName folder and attribute

Kwargs are passed to underlying functions

**Parameters**

- **ExperimentName** (*str*) – Name of experimental ExperimentName
- **Type** (*Optional[str, object]*) – Type of experimental ExperimentName (Optional, default = Experiment)

**Return type**

Any

**create\_log\_file()** → Self

Creates log file

**Return type**

Any

**end\_log()** → Self

Ends Logging

**Return type**

Any

**property experimental\_condition: str**

Experiment condition of the mouse

**Return type**

str

**property instance\_date: str**

Date when this experimental hierarchy was created

**Return type**

str

**classmethod** `load(Directory: str) → Mouse`

Function that loads the entire mouse

**Parameters**

**Directory** (*str*) – Directory containing the organization.json file and associated data

**Returns**

Mouse

**Return type**

ExperimentManagement.Organization.Mouse

**property** `log_file: str`

Log Filename Path

**Return type**

str

**property** `mouse_id: str`

ID of Mouse

**Return type**

str

**pass\_meta()** → Tuple[str, str]

Passes directory/mouse id

**Returns**

directory/mouse id

**Return type**

tuple[str, str]

**record\_experiment\_mod(ExperimentNameKey: str, \*args) → Self**

Record modification of experiment (Data, Time, \*args)

**Parameters**

- **ExperimentNameKey** (*str*) – The key name for the ExperimentName
- **args** (*str*) – A string explaining the modification

**Return type**

Any

**record\_mod(\*args: str) → Self**

Record modification of experiment (Data, Time, \*args)

**Parameters**

**args** (*str*) – A string explaining the modification

**Return type**

Any

**save()** → Self

Saves Mouse to json

**Return type**

Any

**start\_log()** → Self

Starts Log

**Return type**

Any

**update\_all\_folder\_dictionaries()** → Self

This function iterates through all behavioral ExperimentNames to update their folder dictionaries

**Return type**

Any

**class** Organization.Study

Bases: object

Organization.**generate\_read\_me**(*AbsoluteFilePath: str, Text: str*) → None

Generate a read me file

**Parameters**

- **AbsoluteFilePath** (*str*) – Filename path
- **Text** (*str*) – Text inside

**Return type**

None

Organization.**get\_date**()

Organization.**get\_time**()

## 1.2 UserInterfaces module

UserInterfaces.**select\_directory**(*\*\*kwargs*) → str

UserInterfaces.**verbose\_copying**(*src, dst*) → None

## INDICES AND TABLES

- `genindex`
- `modindex`
- `search`





## PYTHON MODULE INDEX

### O

Organization, [1](#)

### U

UserInterfaces, [10](#)

## INDEX

\spxentryadd\_image\_sampling\_folder()\spxextraOrganization.Experiment method, 3  
 \spxentryadd\_notes()\spxextraOrganization.ImagingAnalysis method, 5  
 \spxentryBehavioralExperiment\spxextraOrganization class in Organization, 1  
 \spxentrychannels\spxextraOrganization.Images property, 5  
 \spxentrycheck\_log()\spxextraOrganization.Mouse method, 8  
 \spxentryclean\_up\_compilation()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentryclean\_up\_motion\_correction()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentrycopy\_raw\_imaging\_data()\spxextraOrganization.Experiment method, 4  
 \spxentrycreate()\spxextraOrganization.Mouse method, 8  
 \spxentrycreate\_experiment()\spxextraOrganization.Mouse method, 8  
 \spxentrycreate\_log\_file()\spxextraOrganization.Mouse method, 8  
 \spxentrycurrent\_ExperimentName\spxextraOrganization.ImagingAnalysis property, 6  
 \spxentryData\spxextraOrganization class in Organization, 2  
 \spxentrydefault\_folders()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentryend\_log()\spxextraOrganization.Mouse method, 8  
 \spxentryExperiment\spxextraOrganization class in Organization, 3  
 \spxentryexperimental\_condition\spxextraOrganization.Mouse property, 8  
 \spxentryExperimentName\_id\spxextraOrganization.Experiment property, 3  
 \spxentryexport\_registration\_to\_denoised()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentryFigures\spxextraOrganization class in Organization, 4  
 \spxentryfile\_format\spxextraOrganization.Images property, 5  
 \spxentryfind\_all\_ext()\spxextraOrganization.Data method, 2  
 \spxentryfind\_matching\_files()\spxextraOrganization.Data method, 2  
 \spxentryfolders\spxextraOrganization.Data property, 3  
 \spxentryframes\spxextraOrganization.Images property, 5  
 \spxentrygenerate\_read\_me()\spxextraOrganization module Organization, 10  
 \spxentryget\_date()\spxextraOrganization module Organization, 10  
 \spxentryget\_time()\spxextraOrganization module Organization, 10  
 \spxentryheight\spxextraOrganization.Images property, 5  
 \spxentryImages\spxextraOrganization class in Organization, 5  
 \spxentryimaging\_files\spxextraOrganization.Images property, 5  
 \spxentryImagingAnalysis\spxextraOrganization class in Organization, 5  
 \spxentryinstance\_date\spxextraOrganization.Data property, 3  
 \spxentryinstance\_date\spxextraOrganization.Experiment property, 4  
 \spxentryinstance\_date\spxextraOrganization.Mouse property, 8  
 \spxentryload()\spxextraOrganization.Mouse class method, 8  
 \spxentryload\_cascade\_exports()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentryload\_data()\spxextraOrganization.BehavioralExperiment method, 1  
 \spxentryload\_data()\spxextraOrganization.Experiment method, 4  
 \spxentryload\_fissa\_exports()\spxextraOrganization.ImagingAnalysis method, 6  
 \spxentryload\_image\_unit2()\spxextraOrganization.ImagingAnalysis method, 7  
 \spxentrylog\_file\spxextraOrganization.Mouse property, 9  
 \spxentrymeta\_files\spxextraOrganization.Images property, 5

---

- \spxentrymodule
  - \spxentryOrganization, 1
  - \spxentryUserInterfaces, 10
- \spxentryMouse\spxextraclass in Organization, 7
- \spxentrymouse\_id\spxextraOrganization.Experiment
  - property, 4
- \spxentrymouse\_id\spxextraOrganization.Mouse
  - property, 9
- \spxentrynum\_imaging\_files\spxextraOrganization.Images
  - property, 5
- \spxentrynum\_meta\_files\spxextraOrganization.Images
  - property, 5
- \spxentryOrganization
  - \spxentrymodule, 1
- \spxentrypass\_meta()\spxextraOrganization.Mouse
  - method, 9
- \spxentrypath\spxextraOrganization.Data
  - property, 3
- \spxentryplanes\spxextraOrganization.Images
  - property, 5
- \spxentryrecord\_experiment\_mod()\spxextraOrganization.Mouse
  - method, 9
- \spxentryrecord\_mod()\spxextraOrganization.Experiment
  - method, 4
- \spxentryrecord\_mod()\spxextraOrganization.Mouse
  - method, 9
- \spxentryreindex()\spxextraOrganization.Data
  - method, 3
- \spxentryreorganize\_bruker\_files()\spxextraOrganization.Images
  - method, 5
- \spxentrysave()\spxextraOrganization.Mouse
  - method, 9
- \spxentryselect\_directory()\spxextra module UserInter-
  - faces, 10
- \spxentrystart\_log()\spxextraOrganization.Mouse
  - method, 9
- \spxentryStudy\spxextraclass in Organization, 10
- \spxentryupdate\_all\_folder\_dictionaries()\spxextraOrganization.Mouse
  - method, 10
- \spxentryupdate\_folder\_dictionary()\spxextraOrganization.Experiment
  - method, 4
- \spxentryUserInterfaces
  - \spxentrymodule, 10
- \spxentryverbose\_copying()\spxextra module UserInter-
  - faces, 10
- \spxentryview\_figure()\spxextraOrganization.Figures
  - method, 4
- \spxentrywidth\spxextraOrganization.Images
  - property, 5